

Activity A: Changing Temperature on an Equilibrium Reaction (Starch-Iodine)

A classic test for starch is to add a few drops of iodine solution. A positive test is indicated by a dramatic color change. In this activity, you will investigate the effect of adding or removing temperature as a stress on a reaction at chemical equilibrium. The equilibrium reaction you will investigate is


$$\text{iodine} + \text{starch} \rightleftharpoons \text{starch-iodine complex}$$

Procedure


1. Heat 100 mL water in a 250-mL beaker to 85°C.
2. Place 100 mL water in a second 250-mL beaker and add ice to make a cold water bath.
3. Make a dilute starch solution by filling a pipet with starch solution and placing 5 drops in each well of the reaction tray. Add 10 mL water to each well and stir with a wooden stirrer. Record the color of the starch solution on the reactant side of the equation in step 5 below.

4. Obtain a bottle of iodine-potassium iodide solution from the Materials Station. Add 1 drop of this solution into each well. The middle well will serve as the control well, while wells 1 and 3 will be experimental wells for temperature changes. Measure the temperature of the control well in degrees and record in Data Table 1.

5. Record the color of the starch-iodine complex in the equation below and in the Control well row of Data Table 1.




6. Remove the tops from the two petri dishes and place the bottoms of both dishes under reaction wells 1 and 3 as shown below.



7. To the petri dish under well 3, add 85°C water until it is almost full. Have one partner hold well 3 down into the warm water while the other partner stirs with a thermometer. Note the highest temperature and compare the color to the control tray. Record the temperature and color in Data Table 1. Increasing the temperature causes a stress on this equilibrium. Based upon the colors of the original equilibrium equation above, what is the direction of shift for this new equilibrium? Record this in Data Table 1.
8. Now add ice water to the petri dish under tray 1. Note the lowest temperature and compare the color to that of the control tray. Record both the temperature and color in the data table. Indicate the direction of equilibrium shift for the new reaction.


9. On the basis of your data, include the word "heat" in the correct blank below. (Remember that if you increase heat for a heat-sensitive equilibrium, the reaction shifts to the opposite side to absorb the added heat.)



10. Rinse your three-well reaction tray with water and dry with a paper towel.


Activity B: Changing Concentration on an Equilibrium Reaction (Milk of Magnesia)

Milk of magnesia (MOM) is used as a laxative or an antacid. It is barely soluble in water (only 0.012 g/L). The following equilibrium is established with the small amount that does dissolve:



The stress that you will place on the above equilibrium reaction will be to remove OH^{-} ions by adding

vinegar. The H^{+} ions from the vinegar react with the aqueous OH^{-} ions to form water (This is how milk of magnesia acts as an antacid.)



This reaction can be monitored visually by adding drops of the acid-base indicator, bromthymol blue. The pH range of colors for this indicator is as follows: Yellow < 6.0 ---- green ---- 7.6 > blue. Yellow, then, indicates an acid with pH below 6.0, blue indicates a base with a pH above 7.6, and green indicates a pH between 6.0 and 7.6.

Procedure

1. Add 3 mL of water to wells 1 and 2 in the reaction tray.
2. Stir the milk of magnesia in well 3 of the reagent tray with a wooden stirrer before transferring.
3. Use a pipet to place 15 drops of the MOM suspension in wells 1 and 2 of the reaction tray.
4. Stir each well with the wooden stirrer, get a bromthymol blue dropping bottle from the Materials Station, and add 6 drops to wells 1 and 2.
5. Record the color of both wells and record the pH for this reaction at equilibrium in Data Table 2. Again, the pH range for bromthymol blue is as follows: Yellow > 6.0 ---- green ---- 7.6 > blue. Well 2 will be a control for comparing colors as you add stresses to well 1. Well 3 will remain empty.
6. Fill a pipet with vinegar from well 4 of the reagent tray.
7. Stress 1: Add 10 drops of vinegar to well 1. Stir quickly with the wooden stirrer. Record in Data Table 2 the color change and range of pH for this reaction.
8. Stress 2: Repeat step 7.
9. Stress 3: Repeat step 7.



10. Continue adding additional increments of 10 drops of vinegar and stirring. What eventually happens?
11. In terms of the equation for the solubility of $\text{Mg}(\text{OH})_2$, explain what happens to the equilibrium each time you add vinegar?
12. Explain why there is a limit to how much vinegar you can add until there is no further change in color.

13. Add 10 drops of 1 M NaOH to well 1. Record the color and pH in Data Table 2.
14. Add another 10 drops of 1 M NaOH to well 1. Record the color and pH in Data Table 2.
15. In terms of the $\text{Mg}(\text{OH})_2$ equilibrium equation, explain what happened to the equilibrium as you added 1 M NaOH?
16. Rinse the 4-well reagent tray, the 3-well reaction tray, and the 4 dropping pipets with water.

Activity C: Changing Pressure on an Equilibrium Reaction (Carbonic Acid)

Carbonated water is made by dissolving carbon dioxide gas into water, forming a weak acid, carbonic acid. Carbonic acid then ionizes into H^+ ions and HCO_3^- ions. The equilibrium reactions involved in making

carbonated water are as follows:



The stress that you will manipulate on this reaction is lowering the pressure and observing the shift in the equilibrium as the reaction changes inside a syringe.

Procedure

1. Pour 10 mL of seltzer water into a 50-mL beaker or a small cup.

2. Place 5 drops of Bogen universal indicator into the 10 mL of seltzer water.

3. Take a 10-mL syringe and push the plunger down to 0.0 mL to remove all the air.

4. Place the tip of the syringe into the solution and draw up 3 mL of the solution. Compare the solution color to the Bogen universal indicator chart and record the color and pH of the solution in the Control row of Data Table 3.

5. The pressure on this equilibrium solution is the same as the room atmospheric pressure.

6. Screw the syringe end-cap tightly onto the syringe tip.

7. Pull the plunger up to the 10-mL mark. As you increase the volume, you are decreasing the pressure on your solution. Observe any changes in the solution.

8. Release the plunger. Does the plunger go back to its original starting position of 3 mL?


9. Repeat pulling back to 10 mL and releasing numerous times. Note the changes occurring in the solution.

Compare the color of the solution at the end of several pulls to the Bogen color chart and record color and pH in Data Table 3. Also record the direction of the equilibrium shift that occurs.

10. How does decreasing the pressure affect the initial equilibrium of this reaction?

Activity D: Common Ion Effect

You have a reaction at equilibrium containing 0.5 M copper(II) chloride dissolved in water.



blue

green

You will shift the equilibrium by adding a "common ion" (a compound that has an ion common to the reaction). A blue color indicates a shift left due to an excess of the $\text{Cu}(\text{H}_2\text{O})_6^{2+}$ complex. A green color indicates a shift right due to an excess of the CuCl_4^{2-} complex.

Procedure

Common Ion Effect

1. Place four test tubes in a test tube rack. Use a marker to label the test tubes 1, 2, 3, and 4.

2. Measure 8 mL of copper(II) chloride from the 120-mL bottle. Add 2 mL (40 drops) to each of the

4 test tubes.

3. Keep test tube 1 as a control and reference for color. Record this color in Data Table 4.

4. Weigh 0.3 g of NaCl into one weigh boat, weigh 0.6 g of NaCl into the second boat, and weigh 0.9 g of NaCl into the third boat.

5. Bend one corner of the 0.3 g NaCl boat in half to make a pouring spout. Carefully pour the 0.3 g of NaCl into test tube 2. Gently stir by swirling the test tube.

6. Likewise pour the 0.6 g into test tube 3 and gently swirl to mix.

7. Now pour the 0.9 g into test tube 4 and gently swirl to mix.

8. Compare the colors of test tubes 2, 3, and 4 to that of test tube 1. Record the colors in Data Table 4.

9. Record your answers in Data Table 4

a. What is the "common ion" that you added to test tubes 2, 3, and 4?

Common ion: _____

b. How is the common ion affecting the equilibrium reaction, as the amount of NaCl increases?

Addition of AgNO₃

1. Obtain a bottle of silver nitrate (AgNO₃) from the materials station. Add 20 drops of silver nitrate to test tube 3 and observe what happens. Record the color in Data Table 4. Based upon color, in which direction

does the equilibrium shift? Record in Data Table 4.

2. Write an equation showing what happened when you added the AgNO₃ and explain why the equilibrium


shifted. Record in Data table 4.

Addition of H₂O

1. Use a dropper pipet to add 80–100 drops of water to test tube 4. Gently stir by swirling the test tube. What

happens to the color? Record this color in Data Table 4.

2. Explain this equilibrium shift in terms of the original equilibrium equation:



blue green


3. Keep all four test tube solutions for Activity E.

Activity E: Altering Equilibrium to Achieve Optimum Yield

Assume your lab group is a team of chemical engineers. Your company has a customer that wants to

purchase large quantities of two aqueous complexes that can be produced from the 0.5 M copper(II)

chloride equilibrium reaction investigated in part D.



blue green

The customer wants Solution 1 to have a maximum concentration of Cu(H₂O)₆²⁺(aq) and Solution 2 to have

a maximum concentration of CuCl₄²⁻(aq).

You must design two experimental procedures that would determine the addition or removal of

concentration and temperature stresses for shifting the equilibrium to produce the maximum amount of

each complex in solution. There are no gases in this reaction, so pressure will not be a factor. Have your

teacher approve both procedures. Then, experiment to demonstrate that your procedures will work. Use

control test tube 1 as a reference for color comparisons. Your trial in the test tubes will help your company

decide how to scale the reaction to an industrial level.

Follow your teacher's instructions for proper disposal of all four solutions.